



Handy Cure S™

Finding Good Pain Relief Treatment Therapy was hard till Handy Cure!

Handy Cure S™ is a Class II medical device (certified CE) for pain relief treatments, based on Quantum Medicine (QM) and synergetic effects of 4 energy sources: LLLT (Cold Laser – Class1), Infra-Red, Static electromagnetic field, visible Red Light:

1. **Laser** (an OSRAM unit) at a wavelength of 905 nm drives the photons (light energy) into the target tissue which is absorbed by the haemoglobin, increases the microcirculation, bringing more oxygen to the cells and triggers a faster process of cells renewal. Supports 3 pre-programmed Frequency Modes for different conditions and type of pains (**Variable** between 1000-2500 Hz for All Pain, **50Hz**-for Short-Term Pain, **5Hz** –long term pain).
2. **Infra-Red (IR)** at a wavelength of 875 nm triggers photo physical changes at the cell membrane, improves circulation, increases level of sodium- potassium driven at the cell level hence removing pain signals at the source.
3. **Static magnetic Field** with an induction of 25 -45 mT keeps the photo-activated molecules in the dispersion mode, hence enhancing their energy levels at molecular and cellular levels. It helps in keeping the laser enhanced processes active.
4. **Pulsed Red Light** at 635 nm wavelength – provides a positive therapeutic effect on inflamed tissues.

Handy Cure S™ has 10+ years of field-proven results being used by Chiropractor groups, Osteopaths, Sport Therapists, Physiotherapists, Massage Therapists, Acupuncturists, Rehab/ Elderly Care Centres worldwide for rheumatic & orthopaedic pain relief and sport injuries treatments.

Handy Cure S™ is highly recommended for the following issues:

- Back pain, sciatica, joints, ligaments and tendons pains caused by gout, arthritis, osteoarthritis, coxarthrosis, myositis, other rheumatoid conditions
- Post-traumatic pain, cartilage disorders, radiculitis, ischialgia
- Body joints pain ("Carpal Tunnel" syndrome, "frozen" joints; shoulder, elbow, knee, wrist & heels).
- Spinal diseases, Osteochondritis, Scoliosis
- Sports injuries (Achilles tendonitis, golf & tennis elbow pains, muscles, tendons and ligaments tears and stretching)
- Muscle spasms, neck pain, lower back pain and wound healing

LLLT Laser Class 1M (also called cold laser) is a low-power laser, gentle and safe, with healing properties, the most positive influences on body tissues, so the laser is particularly common used in pain-treatment medicine.

Laser Class1M is more effective than Laser Class 3B in controlling the depth of energy penetration into the tissues and it is safer (no thermal risks, epidermis damage, etc) and it's cheaper (by thousands of dollars less) than Class 4 Lasers (surgical or so called "therapeutic" lasers). Depth of penetration into the tissues body is determined by laser wavelength, NOT power.

LLLT does not produce heat and does not scorch. The healing power of low power laser without any risks or side effects clearly demonstrate a great success. No need for surgery, nor anaesthesia or any drugs.

In addition to pulsating laser energy, **Handy Cure S™** combines synergistically broadband infrared energy, visible red light and static magnetic field, in an unique way that reduces pain and accelerates the formation of new tissue, improving the patients' health, quickly and efficiently and within a relatively short time.

Most painkillers (conventional treatment / topical medication) , mimic the action of endorphins, getting only temporary relief of pain by raising the pain threshold, without actually dealing with the cause of the pain/ inflammation, as actually Handy Cure S™ does.

Healing is de-facto achieved by this device by helping the body's natural mechanisms, namely, improving blood supply to the painful area and improving function at the cellular level (e.g. cartilage cells, bone cells, muscle cells, tendon cells, etc.), improving the blood flow, strengthening the immune system and triggering effective absorption of calcium in cells, improving the blood microcirculation, the cells renewal - these being the true factors for curing the pain/ infection, and for an effective disappearance of the pain.

Referral from one of the US chiro's, one of many [Dr. Craig S. Ross, D.C.,](#)

"I use the Handy Cure in my chiropractic office for many different conditions. I have used more powerful lasers, more expensive lasers -- but the Handy Cure does a fantastic job on my patient's conditions. It is easy to use... so easy that some of my patients have purchased ones for home use. I have 350 doctors in my group and many have purchased the Handy Cure to help even the most difficult patients. Powerful, easy to use and inexpensive is a winning combination".

Handy Cure S™ can be provided with a Set of 4 Optical / Light Acupuncture-Muscle Probes for: Auricular Acupuncture, Muscle Trigger Point Probe, Corporal Probe, Skin/ Wound Scanning Probe (the Optic nozzles are designed to extend the device's capabilities through delivery of a pinpoint laser beam).

